



Luna Auto CPAP System

E-20A-H-0/E-20AJ-H-0

With Luna, you will be fully supported on the path of sleep apnea treatment. E-20A-H-0 and E-20AJ-H-0 assembled with LCD screen and its 3-button control panel simplifies your life.



Luna Auto CPAP System

BMC

E-20A-H-O/E-20AJ-H-O

Easy data access

iCode and SD card are two easy ways to transfer how you went with therapy. Those records can be viewed and evaluated through BMC+ iCode App or PC software.

01

02

Smart pressure match(APAP)

APAP prevents upper airway collapsing with varied therapy pressure. Reslex respiratory pressure relief technology reduces the work of exhale so users remain comfortable and compliance.



Other key features

Ramp | Mask test | Delay off | Multi-language

03

BMC Medical Co., Ltd.

Add: 5/F Main Building, No.19 Gucheng Street West, Shijingshan, Beijing 100043, China

Tel: +86-10-51663880

Fax: +86-10-51663880 ext.810

en.bmc-medical.com

Date: Jun. 12, 2018



User Basic Information

Height(cm):

Weight(kg):

Age:

Sex:

Device Information

Device Serial Number: E2318313238

Model: E-20A

Item	1	7	30	60	90	182	365
Days of Therapy	0	0	0	0	0	0	0
Days of Therapy (>= 4 h)	0	0	0	0	0	0	0
Days of Therapy (%)	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%	0.0%
Avg. Daily Use Time (hh:mm)	00:00	00:00	00:00	00:00	00:00	00:00	00:00
Mean P (cmH ₂ O)	0.0	0.0	0.0	0.0	0.0	0.0	0.0
P95 (cmH ₂ O)	0.0	0.0	0.0	0.0	0.0	0.0	0.0
AHI	0.0	0.0	0.0	0.0	0.0	0.0	0.0
AI	0.0	0.0	0.0	0.0	0.0	0.0	0.0
CAI	0.0	0.0	0.0	0.0	0.0	0.0	0.0
HI	0.0	0.0	0.0	0.0	0.0	0.0	0.0
Avg. Leak (L/min)	0.0	0.0	0.0	0.0	0.0	0.0	0.0
High Leak Time	00:00	00:00	00:00	00:00	00:00	00:00	00:00
SpO ₂ (%)	0.0	0.0	0.0	0.0	0.0	0.0	0.0
ODI	0.0	0.0	0.0	0.0	0.0	0.0	0.0
Res.Rate (BPM)	0.0	0.0	0.0	0.0	0.0	0.0	0.0