

Great Product !

"After reading some of the reviews, I was skeptical that this would be any good. When I took it out of the box, it seemed a little flimsy but was well engineered. Now that I've use it for a week, I wonder why it took so long to decide to buy it. It's a GREAT product that works super. I no longer have a stiff neck from wrestling with the tubing. I LOVE IT!" -LARRY P, TX

## A Real Godsend

"The Hose Lift System was a real godsend. It manages the hose as if it were a real person standing by the bed all night making sure that the hose doesn't interfere with my sleep whenever I turn in bed. I was so delighted with the performance and price of the System, that I ordered another one to have as a spare. I most highly recommend this system to all side sleepers who turn frequently and you can't beat the price." –Leonard R, Arlington, VA

### Don't Hesitate

"I tried 2 hose management systems before this one, and this CPAP Hose Lift is the best -- the easiest to travel with but also the most functional. They key is the flexibility of the pole; the fact that the pole can turn as needed (i.e., if tugged lightly on); and the elastic band at the top, which allows the tube to come down if needed. Really wish I started with this one. -Michael S, NC

Testimonial proven to increase compliant use' by reducing loss of face to mask seal, hose hassle; and, "rain-out" collection with humidifier air. Available from:


Retailer Insert Your Desired Text and Pricing Here

# Ranked #I in Consummer Hose Management Survey



designed of more expensive aluminum

Manufacturer: *Irden Innovations* See: ardeninnovations.com



Hrden Innovations.com



**Suspends your CPAP Tubing Overhead,** Rotates with sleep movement, and thus, Eliminates the Hose Hassle inherent with All CPAP systems.

#### The Essential

# "4th element"

fits within any CPAP machine travel case





'Folding, yet all parts strung gathered in order, it weighs only 5.5 oz (156g) and is a very small 2.1" diameter x 7" long.


# **Dual Use: Home & Travel**

## The new improved CPAP Hose Lift v5.5

is simply the lightest, most compact home and travel, Dual Use hose management device available.

# **Compared to Previous v5**

12% Lighter 16% More Compact 14% Larger Base New: Optional Height / Reach



with v5.5

**Optional height** and reach - just flip the bend, and then the v5.5 will reach out a bit farther and stand a bit shorter. Start by 'unknotting' bungee at arrow:



CE mark: CA011750; HTS#: 9019.20.00; NPI#: 1467684712; Insurance Code#: 1399; Model: CPAPHLv5.5



The first and (only) true Original, with a Curly Cord "roller" on the suspending loop, eases movement of the cpap hose within.



A deep socket mast connection, a mast to base <a href="stretched">stretched</a> cord, "tighten-n-lock:" mast to base - provides an even more stable platform, in hose suspension support.



Includes a Hook for your mask to store, out of your way.